


# FFPE, Platelets & Biospecimens - Cool Cube™ PCM Transport Coolers



A Cool Cube™ at Room Temperatures is not an ordinary transport cooler. It is a portable, validated container and pack-out system that keeps temperature sensitive product (FFPE, platelets, biospecimens, etc.) between 15 and 25 °C for over 3 days (or between 20 and 24 °C for over 2 days) – *without using ice or electricity!* This simple system streamlines the transport of valuable, room temperature product and protects it from warming or cooling during off-site use or power outages. The validated, patented technology of a Cool Cube™ saves time, saves money and is improving cold chain management practices around the world. The Cool Cube™ is a container and pack-out system so advanced...it's easy!



	Cool Cube™ 03 at Room Temps (RT-03)	Cool Cube™ 08 at Room Temps (RT-08)	Cool Cube™ 28 at Room Temps (RT-28)	Cool Cube™ 96 at Room Temps (RT-96)
<b>Temperature Range</b>	20 to 24 °C / 15 to 25 °C	20 to 24 °C / 15 to 25 °C	20 to 24 °C / 15 to 25 °C	20 to 24 °C / 15 to 25 °C
<b>Validated Duration</b>	47+ hrs. / 91+ hrs. (summer)	66+ hrs. / 83+ hrs. (summer)	85+ hrs. / 141+ hrs. (summer)	91+ hrs. / 143+ hrs. (summer)
<b>Payload Area</b>	5¾" x 5¾" x 5¾"	11¼" x 6½" x 6½"	12" x 12" x 12"	18" x 18" x 18"
<b>Exterior Dimensions</b>	11" x 11" x 11"	16" x 11½" x 11½"	20½" x 20½" x 23"	26½" x 26½" x 28¼"
<b>Weight</b>	11 lbs.	17 lbs.	62 lbs.	114 lbs.
<b>Features</b>	<ul style="list-style-type: none"> <li>• Six Room Temp PCM Panels</li> <li>• Temp-Shield™ VIP System</li> <li>• User-Friendly Carrying Case</li> </ul>	<ul style="list-style-type: none"> <li>• Six Room Temp PCM Panels</li> <li>• Temp-Shield™ VIP System</li> <li>• User-Friendly Carrying Case</li> </ul>	<ul style="list-style-type: none"> <li>• Six Room Temp PCM Panels</li> <li>• Temp-Shield™ VIP System</li> <li>• User-Friendly Wheeled Case</li> </ul>	<ul style="list-style-type: none"> <li>• Six Room Temp PCM Panels</li> <li>• Temp-Shield™ VIP System</li> <li>• User-Friendly Wheeled Case</li> </ul>

[Go to Web Page >>](#)



info@VeriCorMed.com  
(608) 526-6901

For more information & pricing: [www.VeriCorMed.com](http://www.VeriCorMed.com)



# Cool Cube™ 03 at Room Temps (FFPE, platelets, biospecimens)

**Model:** RT-03



The Cool Cube™ 03 at Room Temperatures is not an ordinary transport cooler. It is an easy-to-carry, validated container and pack-out system that keeps temperature sensitive product (FFPE, platelets, biospecimens) between 15 and 25 °C for over 3½ days (or between 20 and 24 °C for over 2 days) – without using ice or electricity! This cooler is ideal for transport, but also excels at protecting valuable product during off-site use and power outages, regardless of hot or freezing conditions. The simple, validated, patented technology is the best compact thermal protection product available to maintain the cold chain!

## Features

- 6 patented room PCM panels manage the temp
- Temp-Shield™ vacuum insulated panel (VIP) system
- Rugged, blue carrying case with a shoulder & hand strap
- Validated under ISTA Test Procedure 7D (Int'l Safe Transit Assoc.)

## Specifications

**Temperature Range / Duration**  
20 to 24 °C (68 to 75 °F) / 47+ hrs.  
15 to 25 °C (59 to 77 °F) / 91+ hrs.  
(summer profiles used for hrs.)

**Payload Area**  
5¼" x 5¼" x 5¼"


**Exterior Dimensions & Weight**  
11" x 11" x 11" @ 11 lbs.




- 1 Adjustable shoulder strap with pad is removable.
- 2 Double zipper closure for easy access.
- 3 Clear external pockets for documents and/or Digital Data Logger (DDL).
- 4 Six removable, temperature-controlling, PCM panels ensure a simple, quick & effective pack-out.

## Validated Pack-Out System


**Remove**  
PCM Panels




**Prep**  
PCM Panels



**Assemble**  
PCM Panels



**Pack-Out**  
Product



**Close**  
Bag



[Go to Web Page >>](#)



info@VeriCorMed.com  
(608) 526-6901

For more information & pricing: [www.VeriCorMed.com](http://www.VeriCorMed.com)

**VERICOR** LLC  
Response-Friendly®

# Cool Cube™ 08 at Room Temps (FFPE, platelets, biospecimens)

**Model:** RT-08



The Cool Cube™ 08 at Room Temperatures is not an ordinary transport cooler. It is an easy-to-carry, validated container and pack-out system that keeps temperature sensitive product (FFPE, platelets, biospecimens) between 15 and 25 °C for over 3½ days (or between 20 and 24 °C for over 2½ days) – without using ice or electricity! This cooler is ideal for transport, but also excels at protecting valuable product during off-site use and power outages, regardless of hot or freezing conditions. The simple, validated, patented technology is the best compact thermal protection product available to maintain the cold chain!

## Features

- 6 patented room PCM panels manage the temp
- Temp-Shield™ vacuum insulated panel (VIP) system
- Rugged, blue carrying case with a shoulder & hand strap
- Validated under ISTA Test Procedure 7D (Int'l Safe Transit Assoc.)

## Specifications

**Temperature Range / Duration**  
20 to 24 °C (68 to 75 °F) / 66+ hrs.  
15 to 25 °C (59 to 77 °F) / 83+ hrs.  
(summer profiles used for hrs.)

**Payload Area**  
11¼" x 6½" x 6½"

**Exterior Dimensions & Weight**  
16" x 11½" x 11½" @ 17 lbs.




4 Crêdo™ inside


- 1 Adjustable shoulder strap with pad is removable.
- 2 Double zipper closure for easy access.
- 3 Clear external pockets for documents and/or Digital Data Logger (DDL).
- 4 Six removable, temperature-controlling, PCM panels ensure a simple, quick & effective pack-out.

## Validated Pack-Out System


**Remove**  
PCM Panels




**Prep**  
PCM Panels



**Assemble**  
PCM Panels



**Pack-Out**  
Product



**Close**  
Bag



[Go to Web Page >>](#)



info@VeriCorMed.com  
(608) 526-6901

For more information & pricing: [www.VeriCorMed.com](http://www.VeriCorMed.com)

**VERICOR** LLC  
Response-Friendly®



# Cool Cube™ 28 at Room Temps (FFPE, platelets, biospecimens)

Model: RT-28

The Cool Cube™ 28 at Room Temperatures is not an ordinary transport cooler. It is a lockable, wheeled, qualified container and pack-out system that keeps temperature sensitive product (FFPE, platelets, biospecimens) between 15 and 25 °C for over 5 days (or between 20 and 24 °C for over 3 days) – without using ice or electricity! This cooler is ideal for transport, but also excels at protecting valuable product during off-site use and power outages, regardless of hot or freezing conditions. The simple, validated, patented technology is the best compact thermal protection product available to maintain the cold chain!

## Features

- 6 patented room PCM panels manage the temp
- Temp-Shield™ vacuum insulated panel (VIP) system
- Rugged, military spec case with large lifting handles & wheels
- Validated under ISTA Test Procedure 7D (Int'l Safe Transit Assoc.)

## Specifications

**Temperature Range / Duration**  
20 to 24 °C (68 to 75 °F) / 85+ hrs.  
15 to 25 °C (59 to 77 °F) / 141+ hrs.  
(summer profiles used for hrs.)

**Payload Area**  
12" x 12" x 12"


**Exterior Dimensions & Weight**  
20½" x 20½" x 23" @ 62 lbs.




- 1 Integrated towing handle for easy transport.
- 2 Holder and flexible grommet enables the use of an external temp monitor.
- 3 Tru-Trac System™ - a user-friendly wheel configuration.
- 4 Six removable, temperature-controlling, PCM panels ensure a simple, quick & effective pack-out.

## Validated Pack-Out System


**Remove**  
PCM Panels




**Prep**  
PCM Panels



**Assemble**  
PCM Panels



**Pack-Out**  
Product



**Close**  
Case



[Go to Web Page >>](#)



info@VeriCorMed.com  
(608) 526-6901

For more information & pricing: [www.VeriCorMed.com](http://www.VeriCorMed.com)



# Cool Cube™ 96 at Room Temps (FFPE, platelets, biospecimens)

**Model:** RT-96



The Cool Cube™ 96 at Room Temperatures is not an ordinary transport cooler. It is a lockable, wheeled, qualified container and pack-out system that keeps temperature sensitive product (FFPE, platelets, biospecimens) between 15 and 25 °C for over 5½ days (or between 20 and 24 °C for over 3½ days) – without using ice or electricity! This cooler is ideal for transport, but also excels at protecting valuable product during off-site use and power outages, regardless of hot or freezing conditions. The simple, validated, patented technology is the best compact thermal protection product available to maintain the cold chain!

## Features

- 6 patented room PCM panels manage the temp
- Temp-Shield™ vacuum insulated panel (VIP) system
- Rugged, military spec case with large lifting handles & wheels
- Validated under ISTA Test Procedure 7D (Int'l Safe Transit Assoc.)

## Specifications

**Temperature Range / Duration**  
20 to 24 °C (68 to 75 °F) / 91+ hrs.  
15 to 25 °C (59 to 77 °F) / 143+ hrs.  
(summer profiles used for hrs.)

**Payload Area**  
18" x 18" x 18"


**Exterior Dimensions & Weight**  
26½" x 26½" x 28¼" @ 114 lbs.




- 1 Folding handle for ease of movement.
- 2 Holder and flexible grommet enables the use of an external temp monitor.
- 3 Multi-directional locking wheel system.
- 4 Six removable, temperature-controlling, PCM panels ensure a simple, quick & effective pack-out.

## Validated Pack-Out System


**Remove**  
PCM Panels




**Prep**  
PCM Panels



**Assemble**  
PCM Panels



**Pack-Out**  
Product



**Close**  
Case



[Go to Web Page >>](#)



info@VeriCorMed.com  
(608) 526-6901

For more information & pricing: [www.VeriCorMed.com](http://www.VeriCorMed.com)

