

# Cool Cube™ Vaccine Transporters (Standard & Varicella)

VeriCor's Cool Cube™ Vaccine Transporters boasts superior thermal protection in convenient, rugged, containers ideal for off-site vaccinations. The revolutionary TIC® (Thermal Isolation Chamber) System eliminates the use of ice or electricity when transporting vaccine, yet it breaks down easily for simple storage and preconditioning.



Part #	Temp Range	Validated Hours	Preconditioning Temp Required	Approx Multi-Dose Vials	Weight
VT-50-S	1° - 10°C (34° - 50°F)	96+	-18°C (0°F)	50	10 lbs
VT-50-V	< -15°C (5°F)	72+	-25°C (-13°F)	50	10 lbs
VT-250-S	1° - 10°C (34° - 50°F)	96+	-18°C (0°F)	250	47 lbs
VT-250-V	< -15°C (5°F)	96+	-25°C (-13°F)	250	47 lbs
VT-420-S	1° - 10°C (34° - 50°F)	96+	-18°C (0°F)	420	58 lbs
VT-420-V	< -15°C (5°F)	96+	-25°C (-13°F)	420	58 lbs



Vaccine Transporters  
TIC® System

