

# Westcot Surge Beds

The Westcot line of surge beds is reliable, safe, and comfortable. This “gold standard” of surge beds are what others try to copy. Beds are tested 10% above rated capacity for a minimum of 5,000 cycles. There are no holes, cracks or sharp edges to cut fingers or allow contaminants to lodge. Meanwhile the superior comfort approaches box spring construction.

## Features

- Life cycle tested
- Aluminum folding frame
- Spun tube ends
- Anti-microbial mattress cover
- Springs connect deck to frame
- Meets federal flammability standard
- Response-Friendly™

## Bed Line

Westcot 400  
Westcot 400 APC  
MedBed 400  
Bariatric Bed 700  
Bio-Hazard De-Con Cot  
Commander Cot  
Surge Cot

## Accessories

Racks  
I.V. Poles  
Bed Safety Rails  
Hanging Drawers  
Rack Dust Covers  
Disposable Linens  
*+ more beds/accessories available*



## Westcot Surge Beds

Part #: call for models & pricing



**VERICOR** LLC  
Response-Friendly™ Medical Systems