

VeriCor's 4-Drawer Medical Case is a storage container, cabinet and mobile case all-in-one. This durable case provides deep widths and easy access to vital medical gear in the field. It is designed specifically to keep life saving equipment/supplies safe and secure from moisture, dust, sand, temperature and impact.


Features

- Wheeled for easy mobility
- Spring-loaded butterfly latches
- Drawers are removable and interchangeable
- Rotationally-molded
- Airtight, waterproof, and sand proof
- Chemical agent decontaminable
- Stackable
- Response-Friendly™



Specifications

Weight: 49 lbs.
Dimensions: 34" x 21" x 19"
Drawer Size (each): 16" x 14.75" x 6.5"



4-Drawer Medical Case
Part #: MC-4



VERICOR LLC
Response-Friendly™ Medical Systems