

VeriCor's 42" Vaccination Workstation comes standard with ample storage, a built-in refrigerator and a detachable table. This single unit station becomes mobile in an instant. Pack in all your medical supplies and the workstation is ready for your next vaccination campaign.



*Features*

- Drawers are constructed with stainless steel for easy decontamination and high durability (partitions included)
- Medical grade drawer pulls
- Locking lid cover stays for safe operations
- Lockable case and refrigerator
- Refrigerator maintains consistent temperature (does not have freezer)
- Built-in power strip with easy access exterior connection
- Heavy duty tie-downs for secure transport
- Wide wheel base
- Available in: Blue, Black, Red, Silver
- Response-Friendly™

*\*complimentary labeling included*



*Specifications*

Dimensions (l x w x h): 43" x 24" x 48"  
 Weight: 270 lbs.  
 Drawers: 3 (100 lbs. capacity each)  
 Fold-Out Table: 8 ft.<sup>2</sup> work surface  
 Tilt Bins: 12 (each 9.5" x 3.5")

*Refrigerator Specifications*

Dimensions (l x w x h): 19" x 18" x 32.5"  
 Capacity: 3.1 ft<sup>3</sup>  
 Voltage: 110V/60Hz\*

*\*220V/50Hz available*



**42" Vaccination Workstation**  
 Part #: MC-42-3-RF-A



**VERICOR** LLC  
 Response-Friendly™ Medical Systems