

# 42" 6-Drawer Workstation

# MC-42-6

VeriCor's 42" 6-Drawer Workstation comes standard with plenty of storage and a detachable table for an instant work surface. This single unit workstation becomes mobile in an instant. Pack in all your medical supplies and the workstation is ready for your next campaign.



### Features

- Drawers are constructed with stainless steel for easy decontamination and high durability (partitions included)
- Medical grade drawer pulls
- Locking lid cover stays for safe operations
- Lockable case for secure storage
- Heavy duty tie-downs for secure transport
- Wide wheel base
- Available in: Blue, Black, Red, Silver
- Response-Friendly™

*\*complimentary labeling included*



### Specifications


Dimensions (l x w x h): 43" x 24" x 48"

Weight: 250 lbs.

Drawers: 6 (100 lbs. capacity each)

Fold-Out Table: 8 ft.<sup>2</sup> work surface

Tilt Bins: 12 (each 9.5" x 3.5" x 5")



## 42" Vaccination Workstation

Part #: MC-42-6



# VERICOR<sup>LLC</sup>

Response-Friendly™ Medical Systems