

48" Vaccination Workstation

MC-48-3-RF

VeriCor's 48" Vaccination Workstation comes standard with ample storage, a built-in refrigerator and a detachable table. This single unit station becomes mobile in an instant. Pack in all your medical supplies and the workstation is ready for your next vaccination campaign.



Features

- Drawers are constructed with stainless steel for easy decontamination and high durability (partitions included)
- Medical grade drawer pulls
- Locking lid cover stays for safe operations
- Lockable case and refrigerator
- Refrigerator maintains consistent temperature (does not have freezer)
- Built-in power strip with easy access exterior connection
- Heavy duty tie-downs for secure transport
- Wide wheel base
- Available in: Blue, Black, Red, Silver
- Response-Friendly™

**complimentary labeling included*



Specifications

Dimensions (l x w x h): 48" x 24" x 48"

Weight: 300 lbs.

Drawers: 3 (100 lbs. capacity each)

Fold-Out Table: 10.4 ft.² work surface

Tilt Bins: 12 (each 9.5" x 3.5")

Refrigerator Specifications

Dimensions (l x w x h): 19" x 18" x 32.5"

Capacity: 3.1 ft³

Voltage: 110V/60Hz*

**220V/50Hz available*



48" Vaccination Workstation

Part #: MC-48-3-RF



VERICOR^{LLC}

Response-Friendly™ Medical Systems