

Certified Traceable Thermometer

VeriCor's Certified Traceable Thermometer accurately monitors temperature in Vaccine Transporters and refrigerators with the attached temperature-buffered sensor. This patented unit fulfills all CDC thermometer and vaccine thermometer requirements.

Features

- Records high/low readings
- Programmable alarm signals out-of-range conditions (visual & audible)
- Avoids contaminants (probe is sealed)
- Traceable to NIST for accuracy
- Response-Friendly™

Specifications

Temperature Range: -50° to 70°C (-58° to 158°F)

Weight: 4 oz.

Included Contents: Certificate, Flip-Open Stand, Wall Mount Slot, 10' Micro-Cable, Velcro®, and Magnetic Strips

V-TCT



Certified Traceable Thermometer

Part #: V-TCT



VERICOR^{LLC}

Response-Friendly™ Medical Systems

info@VeriCorMed.com

(866) 469-6019

www.VeriCorMed.com