

The Westcart Bed Rolling Rack protects Westcot beds during storage and deployment. It moves easily on carpet, vinyl, wood and concrete flooring. It features a welded steel construction with laser cut components for no sharp edges.



Features

- Welded steel construction
- Laser cut components
- Baked-on, powder-coated, textured finish
- Two fixed casters, two swivel casters
- 8 or 10-bed configurations*
- Response-Friendly™

**8-bed configuration shown; does not include beds*

Specifications

Frame Material:	Laser cut steel
Rack Footprint:	38" x 43"
8 Bed Rack Height:	70" (for trailers)
10 Bed Rack Height:	82"
Weight:	155 lbs.
Casters:	3" diameter

shown: 3 Racks/25 Beds, standard in the Medical Equipment Module



Westcart Bed Rolling Rack
Part #: WCR-8 or WCR-10



VERICOR LLC
Response-Friendly™ Medical Systems