

5L Medical Equipment Response Case

Model: MC-MERC-005L

The 5L Medical Equipment Response Case provides maximum protection for contents, making it ideal for valuable medical equipment. The airtight, unbreakable, Pelican case, combined with a tamper seal lock and custom labeling, makes storing and transporting equipment quick and easy. Pack up the equipment you value most into this Response-Friendly™ case and be ready to respond.

Customizable Foam Interior



Lock Included



Sample Labeling



Features

- Strong, lightweight, solid walls
- Pick N Pluck™ foam
- Easy-open, double-throw latches
- Complimentary alphanumeric identification labeling
- Lifetime guarantee

Specifications

Dimensions (L x W x D)
11" x 10" x 5"

Interior (L x W x D)
9¼" x 7" x 4"

Weight
3 lbs.

Interior Cubic Volume
272 cubic inches

Used in Response-Friendly™ modules.

[Go to Web Page >>](#)



SS-MCMERC005L-022719

info@VeriCorMed.com
(866) 469-6019

For more information & pricing: www.VeriCorMed.com

VERICOR LLC
Response-Friendly®