

# 32L Medical Equipment Response Case

**Model:** MC-MERC-032L

The 32L Medical Equipment Response Case provides maximum protection for contents, making it ideal for valuable medical equipment. The airtight, unbreakable, Pelican case, combined with a tamper seal lock and custom labeling, makes storing and transporting equipment quick and easy. Pack up the equipment you value most into this Response-Friendly™ case and be ready to respond.

Customizable Foam Interior



Lock Included



## Features

- Strong, lightweight, solid walls
- Pick N Pluck™ foam
- Easy-open, double-throw latches
- Complimentary alphanumeric identification labeling
- Lifetime guarantee

## Specifications

*Dimensions (L x W x D)*  
21" x 18" x 9"

*Interior (L x W x D)*  
18¾" x 14" x 7¾"

*Weight*  
12 lbs.

*Interior Cubic Volume*  
2040 cubic inches

*Used in Response-Friendly™ modules.*



Sample Labeling



[Go to Web Page >>](#)



SS-MCMERC032L-022719

info@VeriCorMed.com  
(866) 469-6019

For more information & pricing: [www.VeriCorMed.com](http://www.VeriCorMed.com)

 **VERICOR** LLC  
Response-Friendly®