

# Fresh Vaccine & Medicine - Cool Cube™ Transport Coolers



A Cool Cube™ at Refrigerator Temperatures is not an ordinary transport cooler. It is a portable, qualified container and pack-out system that keeps temperature sensitive product (fresh vaccine, medicine, etc.) chilled between 2° and 8°C for over 2.5 days – *without using ice or electricity!* This simple system streamlines the transport of valuable, refrigerated product and protects it from warming or freezing during off-site use or power outages. The validated, patented technology of a Cool Cube™ saves time, saves money and is improving cold chain management practices around the world. The Cool Cube™ is a container and pack-out system so advanced...it's easy!



	Cool Cube™ 03 at Refrigerator Temps (VT-03)	Cool Cube™ 08 at Refrigerator Temps (VT-08)	Cool Cube™ 28 at Refrigerator Temps (VT-28)	Cool Cube™ 96 at Refrigerator Temps (VT-96)
<b>Temperature Parameter</b>	2° to 8°C (36° to 46°F)	2° to 8°C (36° to 46°F)	2° to 8°C (36° to 46°F)	2° to 8°C (36° to 46°F)
<b>Validated Duration</b>	65+ hrs.	76+ hrs.	103+ hrs.	126+ hrs.
<b>Payload Dimensions</b>	5¾" x 5¾" x 5¾"	16" x 11½" x 11½"	12" x 12" x 12"	18" x 18" x 18"
<b>Exterior Dimensions</b>	11" x 11" x 11"	11¼" x 6½" x 6½"	20½" x 20½" x 23"	26½" x 26½" x 28¼"
<b>Weight</b>	11 lbs.	16 lbs.	61 lbs.	111 lbs.
<b>Features</b>	<ul style="list-style-type: none"> <li>• 6 PCM Panels</li> <li>• Temp-Shield™ VIP System</li> <li>• User-Friendly Carrying Case</li> </ul>	<ul style="list-style-type: none"> <li>• 6 PCM Panels</li> <li>• Temp-Shield™ VIP System</li> <li>• User-Friendly Carrying Case</li> </ul>	<ul style="list-style-type: none"> <li>• 6 PCM Panels</li> <li>• Temp-Shield™ VIP System</li> <li>• User-Friendly Wheeled Case</li> </ul>	<ul style="list-style-type: none"> <li>• 6 PCM Panels</li> <li>• Temp-Shield™ VIP System</li> <li>• User-Friendly Wheeled Case</li> </ul>

[Go to Web Page >>](#)



info@VeriCorMed.com  
(866) 469-6019

For more information & pricing: [www.VeriCorMed.com](http://www.VeriCorMed.com)

