


# FFPE, Platelets & Biospecimens - Cool Cube™ Transport Coolers



A Cool Cube™ at Controlled Room Temperatures is not an ordinary transport cooler. It is a portable, qualified container and pack-out system that keeps temperature sensitive product (FFPE, platelets, biospecimens, etc.) between 15° and 25°C for over 3 days – *without using ice or electricity!* This simple system streamlines the transport of valuable, room temperature product and protects it from warming or cooling during off-site use or power outages. The validated, patented technology of a Cool Cube™ saves time, saves money and is improving cold chain management practices around the world. The Cool Cube™ is a container and pack-out system so advanced...it's easy!



	Cool Cube™ 03 at Controlled Room Temps (RT-03)	Cool Cube™ 08 at Controlled Room Temps (RT-08)	Cool Cube™ 28 at Controlled Room Temps (RT-28)	Cool Cube™ 96 at Controlled Room Temps (RT-96)
<b>Temperature Parameter</b>	15° to 25°C (59° to 77°F)	15° to 25°C (59° to 77°F)	15° to 25°C (59° to 77°F)	15° to 25°C (59° to 77°F)
<b>Validated Duration</b>	47+ hrs. (summer)	66+ hrs. (summer)	79+ hrs. (summer)	132+ hrs. (summer)
<b>Payload Dimensions</b>	5¾" x 5¾" x 5¾"	11¼" x 6½" x 6½"	12" x 12" x 12"	18" x 18" x 18"
<b>Exterior Dimensions</b>	11" x 11" x 11"	16" x 11½" x 11½"	20½" x 20½" x 23"	26½" x 26½" x 28¼"
<b>Weight</b>	11 lbs.	17 lbs.	62 lbs.	114 lbs.
<b>Features</b>	<ul style="list-style-type: none"> <li>• 6 PCM Panels (Crêdo™ Series 22)</li> <li>• Temp-Shield™ VIP System</li> <li>• User-Friendly Carrying Case</li> </ul>	<ul style="list-style-type: none"> <li>• 6 PCM Panels (Crêdo™ Series 22)</li> <li>• Temp-Shield™ VIP System</li> <li>• User-Friendly Carrying Case</li> </ul>	<ul style="list-style-type: none"> <li>• 6 PCM Panels (Crêdo™ Series 22)</li> <li>• Temp-Shield™ VIP System</li> <li>• User-Friendly Wheeled Case</li> </ul>	<ul style="list-style-type: none"> <li>• 6 PCM Panels (Crêdo™ Series 22)</li> <li>• Temp-Shield™ VIP System</li> <li>• User-Friendly Wheeled Case</li> </ul>

[Go to Web Page >>](#)



info@VeriCorMed.com  
(866) 469-6019

For more information & pricing: [www.VeriCorMed.com](http://www.VeriCorMed.com)

