

Varicella, MMRV, Zoster & FFP - Cool Cube™ Transport Coolers

A Cool Cube™ at Frozen Temperatures is not an ordinary transport cooler. It is a portable, qualified container and pack-out system that keeps temperature sensitive product (varicella, MMRV, zoster, Fresh Frozen Plasma) frozen between -50° and -15°C for over 2.5 days – without using ice or electricity!

Temperature Range | -50° to -15°C (-58° to 5°F)



	Cool Cube™ 03 at Frozen Temps (FT-03)	Cool Cube™ 08 at Frozen Temps (FT-08)	Cool Cube™ 28 at Frozen Temps (FT-28)	Cool Cube™ 96 at Frozen Temps (FT-96)
Validated Duration	62+ hrs.	60+ hrs.	94+ hrs.	139+ hrs.
Payload Dimensions	5.75" x 5.75" x 5.75"	11.25" x 6.5" x 6.5"	12" x 12" x 12"	18" x 18" x 18"
Exterior Dimensions	11" x 11" x 11"	16" x 11.5" x 11.5"	21" x 22" x 23"	27" x 26" x 25"
Weight	13 lbs.	21 lbs.	71 lbs.	123 lbs.
Features	<ul style="list-style-type: none"> • 6 PCM Panels (Crêdo™ Series 20M) • High-Performance VIP System • User-Friendly Carrying Case 	<ul style="list-style-type: none"> • 6 PCM Panels (Crêdo™ Series 20M) • High-Performance VIP System • User-Friendly Carrying Case 	<ul style="list-style-type: none"> • 6 PCM Panels (Crêdo™ Series 20M) • High-Performance VIP System • User-Friendly Wheeled Case 	<ul style="list-style-type: none"> • 6 PCM Panels (Crêdo™ Series 20M) • High-Performance VIP System • User-Friendly Wheeled Case
	Web/Pricing	Web/Pricing	Web/Pricing	Web/Pricing