


Blood Products - Cool Cube™ Transport Coolers



A Cool Cube™ at Fridge Temperatures is not an ordinary transport cooler. It is a portable, qualified container and pack-out system that keeps temperature sensitive product (whole blood, RBCs, liquid plasma, etc.) chilled between 1° and 10°C for over 2.5 days – *without using ice or electricity!* This simple system streamlines the transport of valuable, refrigerated product and protects it from warming or freezing during off-site use or power outages. The validated, patented technology of a Cool Cube™ saves time, saves money and is improving cold chain management practices around the world. The Cool Cube™ is a container and pack-out system so advanced...it's easy!



	Cool Cube™ 03 at Fridge Temps (BT-03)	Cool Cube™ 08 at Fridge Temps (BT-08)	Cool Cube™ 28 at Fridge Temps (BT-28)	Cool Cube™ 96 at Fridge Temps (BT-96)
Temperature Parameter	1° to 10°C (34° to 50°F)	1° to 10°C (34° to 50°F)	1° to 10°C (34° to 50°F)	1° to 10°C (34° to 50°F)
Validated Duration	70+ hrs.	83+ hrs.	108+ hrs.	128+ hrs.
Payload Dimensions	5¾" x 5¾" x 5¾"	16" x 11½" x 11½"	12" x 12" x 12"	18" x 18" x 18"
Exterior Dimensions	11" x 11" x 11"	11¼" x 6½" x 6½"	20½" x 20½" x 23"	26½" x 26½" x 28¼"
Weight	11 lbs.	16 lbs.	61 lbs.	111 lbs.
Features	<ul style="list-style-type: none"> • 6 PCM Panels (Crēdo™ Series 4) • Temp-Shield™ VIP System • User-Friendly Carrying Case 	<ul style="list-style-type: none"> • 6 PCM Panels (Crēdo™ Series 4) • Temp-Shield™ VIP System • User-Friendly Carrying Case 	<ul style="list-style-type: none"> • 6 PCM Panels (Crēdo™ Series 4) • Temp-Shield™ VIP System • User-Friendly Wheeled Case 	<ul style="list-style-type: none"> • 6 PCM Panels (Crēdo™ Series 4) • Temp-Shield™ VIP System • User-Friendly Wheeled Case

[Go to Web Page >>](#)



info@VeriCorMed.com
(866) 469-6019

For more information & pricing: www.VeriCorMed.com

