

48" Workstation w/ 3 Drawers, Table & Refrigerator

Model: MW-48R

VeriCor's 48" Workstation with Refrigerator, 3 Drawers & Table comes standard with ample storage, a built-in refrigerator and a detachable table. This single unit station becomes mobile in an instant. Pack in all your medical supplies and the workstation is ready for your next vaccination campaign.

Integrated Refrigerator



Locking Lid Stays



Swivel Wheels



Features

- Refrigerator without freezer
- Lockable refrigerator & case
- In-lid glove holders
- Easy lock wheels with a wide wheel base for easy and safe maneuverability
- Custom labeling included

Specifications

Dimensions (H x W x D)
48" x 49" x 25"

Drawers
3 (100 lbs. capacity each)

Fold-Out Table
35" x 49" (30" tall)

Weight
323 lbs.

Please provide color code when ordering:

COLOR (CODE)	Red (RD)	Yellow (YL)	Forest Green (FG)	Navy Blue (NB)	Black (BK)
Silver (SL)	Orange (OR)	Green (GR)	Violet Blue (BL)	Wine (WN)	Tan (TN)

Go to Web Page >>



info@VeriCorMed.com
(866) 469-6019

For more information & pricing: www.VeriCorMed.com

VERICOR LLC
Response-Friendly®