


Active Patient Care Bed w/ IV Pole

Model: RB-APC400I

The Active Patient Care (APC) Bed is designed with patient surge treatment care needs in mind. This taller, procedure height bed allows healthcare professionals to perform necessary functions at a working height. In the event of an overflow, the APC doubles effectively as a temporary cot.



Folded (Dual View)



Features

- Lightweight and portable
- 5 position head elevator
- 2 position foot elevator
- Anti-microbial mattress
- Spring and deck design
- Patient restraint belt
- Stow on board IV Pole

Specifications

- 450+ lbs. working capacity
- Sleeping Surface: 32" x 80"
- Height: 28"
- Folded: 34" x 42" x 10"
- Weight: 38 lbs.
- Meets FED STD 1633 for flammability

| Options | Model # |
|--------------------|--------------|
| w/ IV Pole & Rails | RB-APC400ISR |

[Go to Web Page >>](#)



SS-RB-APC400I-012919

info@VeriCorMed.com
(866) 469-6019

For more information & pricing: www.VeriCorMed.com

VERICOR LLC
Response-Friendly®