

# Bariatric Bed, Wt Cap 700lb

**Model:** RB-B700

The Bariatric Response Bed is a lightweight, portable, reusable cot that is highly versatile. The extra sturdy design offers the bariatric population a safe and comfortable bed in times of need. In addition, the taller, longer design has been found to accommodate multiple users, maximizing the limited space of a patient surge situation.



Folded (Dual View)



## Features

- Lightweight and portable
- Spun tube ends
- Waterproof, anti-microbial mattress
- Comfortable spring and deck design
- Heavy duty locking legs

## Specifications

- 700+ lbs. working capacity
- Sleeping Surface: 38" x 84"
- Height: 22"
- Folded: 42" x 42" x 7"
- Weight: 36 lbs.
- Meets FED STD 1633 for flammability



*\* linen set sold separately*

[Go to Web Page >>](#)



info@VeriCorMed.com  
(866) 469-6019

For more information & pricing: [www.VeriCorMed.com](http://www.VeriCorMed.com)

**VERICOR** LLC  
Response-Friendly®