


Cool Cube™ PCM Coolers

Easy, Qualified, Pack-Out System w/ Passive Cooling Technology



5 Simple Steps



@ Refrigerator Temps
for Fresh Vaccine & Medicine

2 to 8 °C

[More Info](#)

@ Fridge Temps
for Blood Products

1 to 6 °C
1 to 10 °C

[More Info](#)

@ Room Temps
for FFPE, Platelets & Biospecimens

15 to 25 °C
20 to 24 °C

[More Info](#)

@ Lab Freezer Temps
for Varicella, MMRV, Zoster & FFP

-50 to -15 °C

[More Info](#)

Cool Cube™ PCM Coolers for Fresh Vaccine & Medicine

A Cool Cube™ at Refrigerator Temperatures is not an ordinary transport cooler. It is a portable, qualified container and pack-out system that keeps temperature sensitive product (fresh vaccine, medicine, etc.) chilled between 2 and 8 °C for over 2½ days – *without using ice or electricity!*

Temperature Range | 2 to 8 °C (36 to 46 °F)



	Cool Cube™ 03 at Refrigerator Temps (VT-03)	Cool Cube™ 08 at Refrigerator Temps (VT-08)	Cool Cube™ 28 at Refrigerator Temps (VT-28)	Cool Cube™ 96 at Refrigerator Temps (VT-96)
Qualified Duration	65+ hrs.	76+ hrs.	103+ hrs.	126+ hrs.
Payload Area	5¼" x 5¼" x 5¼"	11¼" x 6½" x 6½"	12" x 12" x 12"	18" x 18" x 18"
Exterior Dimensions	11" x 11" x 11"	16" x 11½" x 11½"	20½" x 20½" x 23"	26½" x 26½" x 28¼"
Weight	11 lbs.	16 lbs.	61 lbs.	111 lbs.
Features	<ul style="list-style-type: none"> • Six Refrigerator PCM Panels • Temp-Shield™ VIP System • User-Friendly Carrying Case 	<ul style="list-style-type: none"> • Six Refrigerator PCM Panels • Temp-Shield™ VIP System • User-Friendly Carrying Case 	<ul style="list-style-type: none"> • Six Refrigerator PCM Panels • Temp-Shield™ VIP System • User-Friendly Wheeled Case 	<ul style="list-style-type: none"> • Six Refrigerator PCM Panels • Temp-Shield™ VIP System • User-Friendly Wheeled Case
	Web/Pricing	Web/Pricing	Web/Pricing	Web/Pricing

Cool Cube™ PCM Coolers for Blood Products

A Cool Cube™ at Fridge Temperatures is not an ordinary transport cooler. It is a portable, validated container and pack-out system that keeps temperature sensitive product (whole blood, RBCs, liquid plasma, etc.) chilled between 1 and 10 °C for over 2½ days (or between 1 and 6 °C for over 1½ days) – *without using ice or electricity!*

Temperature Range | 1 to 6 °C (34 to 42 °F) / 1 to 10 °C (34 to 50 °F)



	Cool Cube™ 03 at Fridge Temps (BT-03)	Cool Cube™ 08 at Fridge Temps (BT-08)	Cool Cube™ 28 at Fridge Temps (BT-28)	Cool Cube™ 96 at Fridge Temps (BT-96)
Validated Duration	39+ hrs. / 70+ hrs.	53+ hrs. / 83+ hrs.	68+ hrs. / 108+ hrs.	112+ hrs. / 128+ hrs.
Payload Area	5¼" x 5¼" x 5¼"	11¼" x 6½" x 6½"	12" x 12" x 12"	18" x 18" x 18"
Exterior Dimensions	11" x 11" x 11"	16" x 11½" x 11½"	20½" x 20½" x 23"	26½" x 26½" x 28¼"
Weight	11 lbs.	16 lbs.	61 lbs.	111 lbs.
Features	<ul style="list-style-type: none"> • Six Refrigerator PCM Panels • Temp-Shield™ VIP System • User-Friendly Carrying Case 	<ul style="list-style-type: none"> • Six Refrigerator PCM Panels • Temp-Shield™ VIP System • User-Friendly Carrying Case 	<ul style="list-style-type: none"> • Six Refrigerator PCM Panels • Temp-Shield™ VIP System • User-Friendly Wheeled Case 	<ul style="list-style-type: none"> • Six Refrigerator PCM Panels • Temp-Shield™ VIP System • User-Friendly Wheeled Case
	Web/Pricing	Web/Pricing	Web/Pricing	Web/Pricing

Cool Cube™ PCM Coolers for FFPE, Platelets & Biospecimens

A Cool Cube™ at Room Temperatures is not an ordinary transport cooler. It is a portable, validated container and pack-out system that keeps temperature sensitive product (FFPE, platelets, biospecimens, etc.) between 15 and 25 °C for over 3 days (or between 20 and 24 °C for over 2 days) – *without using ice or electricity!*

Temperature Range | 20 to 24 °C (68 to 75 °F) / 15 to 25 °C (59 to 77 °F)



	Cool Cube™ 03 at Room Temps (RT-03)	Cool Cube™ 08 at Room Temps (RT-08)	Cool Cube™ 28 at Room Temps (RT-28)	Cool Cube™ 96 at Room Temps (RT-96)
Validated Duration	47+ hrs. / 91+ hrs. (summer)	66+ hrs. / 83+ hrs. (summer)	85+ hrs. / 141+ hrs. (summer)	91+ hrs. / 143+ hrs. (summer)
Payload Area	5¼" x 5¼" x 5¼"	11¼" x 6½" x 6½"	12" x 12" x 12"	18" x 18" x 18"
Exterior Dimensions	11" x 11" x 11"	16" x 11½" x 11½"	20½" x 20½" x 23"	26½" x 26½" x 28¼"
Weight	11 lbs.	17 lbs.	62 lbs.	114 lbs.
Features	<ul style="list-style-type: none"> • Six Room Temp PCM Panels • Temp-Shield™ VIP System • User-Friendly Carrying Case 	<ul style="list-style-type: none"> • Six Room Temp PCM Panels • Temp-Shield™ VIP System • User-Friendly Carrying Case 	<ul style="list-style-type: none"> • Six Room Temp PCM Panels • Temp-Shield™ VIP System • User-Friendly Wheeled Case 	<ul style="list-style-type: none"> • Six Room Temp PCM Panels • Temp-Shield™ VIP System • User-Friendly Wheeled Case
	Web/Pricing	Web/Pricing	Web/Pricing	Web/Pricing

Cool Cube™ PCM Coolers for Varicella, MMRV, Zoster & FFP

A Cool Cube™ at Lab Freezer Temperatures is not an ordinary transport cooler. It is a portable, qualified container and pack-out system that keeps temperature sensitive product (varicella, MMRV, zoster, Fresh Frozen Plasma) frozen between -50 and -15 °C for over 2½ days – without using ice or electricity!

Temperature Range | -50 to -15 °C (-58 to 5 °F)



	Cool Cube™ 03 at Lab Freezer Temps (FT-03)	Cool Cube™ 08 at Lab Freezer Temps (FT-08)	Cool Cube™ 28 at Lab Freezer Temps (FT-28)	Cool Cube™ 96 at Lab Freezer Temps (FT-96)
Qualified Duration	62+ hrs.	60+ hrs.	94+ hrs.	139+ hrs.
Payload Area	5¼" x 5¼" x 5¼"	11¼" x 6½" x 6½"	12" x 12" x 12"	18" x 18" x 18"
Exterior Dimensions	11" x 11" x 11"	16" x 11½" x 11½"	20½" x 20½" x 23"	26½" x 26½" x 28¼"
Weight	14 lbs.	20 lbs.	70 lbs.	134 lbs.
Features	<ul style="list-style-type: none"> • Six Lab Freezer PCM Panels • Temp-Shield™ VIP System • User-Friendly Carrying Case 	<ul style="list-style-type: none"> • Six Lab Freezer PCM Panels • Temp-Shield™ VIP System • User-Friendly Carrying Case 	<ul style="list-style-type: none"> • Six Lab Freezer PCM Panels • Temp-Shield™ VIP System • User-Friendly Wheeled Case 	<ul style="list-style-type: none"> • Six Lab Freezer PCM Panels • Temp-Shield™ VIP System • User-Friendly Wheeled Case
	Web/Pricing	Web/Pricing	Web/Pricing	Web/Pricing