

FFPE, Platelets & Biospecimens - Cool Cube™ Transport Coolers

A Cool Cube™ at Controlled Room Temperatures is not an ordinary transport cooler. It is a portable, qualified container and pack-out system that keeps temperature sensitive product (FFPE, platelets, biospecimens, etc.) between 15° and 25°C for over 3 days – *without using ice or electricity!*

Temperature Range | 15° to 25°C (59° to 77°F) / 20° to 24°C (68° to 75°F)



	Cool Cube™ 03 at Controlled Room Temps (RT-03)	Cool Cube™ 08 at Controlled Room Temps (RT-08)	Cool Cube™ 28 at Controlled Room Temps (RT-28)	Cool Cube™ 96 at Controlled Room Temps (RT-96)
Validated Duration	91+ hrs. (summer)	83+ hrs. (summer)	79+ hrs. (summer)	132+ hrs. (summer)
Payload Dimensions	5.75" x 5.75" x 5.75"	11.25" x 6.5" x 6.5"	12" x 12" x 12"	18" x 18" x 18"
Exterior Dimensions	11" x 11" x 11"	16" x 11.5" x 11.5"	21" x 22" x 23"	27" x 26" x 25"
Weight	11 lbs.	17 lbs.	62 lbs.	105 lbs.
Features	<ul style="list-style-type: none"> • 6 PCM Panels (Crêdo™ Series 22) • High-Performance VIP System • User-Friendly Carrying Case 	<ul style="list-style-type: none"> • 6 PCM Panels (Crêdo™ Series 22) • High-Performance VIP System • User-Friendly Carrying Case 	<ul style="list-style-type: none"> • 6 PCM Panels (Crêdo™ Series 22) • High-Performance VIP System • User-Friendly Wheeled Case 	<ul style="list-style-type: none"> • 6 PCM Panels (Crêdo™ Series 22) • High-Performance VIP System • User-Friendly Wheeled Case
	Web/Pricing	Web/Pricing	Web/Pricing	Web/Pricing