

# Emergency Room Bed

**Model:** RB-ER

The Emergency Room Response Bed is the ideal bed for mobile emergency rooms. The heavy-duty construction withstands the rigors of field medicine, while the features enable safe, efficient treatment. Push-button assembly allows for a quick setup without tools. After use, all parts can be easily disinfected and decontaminated for future use.



## Features

- Hard deck for rigorous medical procedures
- Adjustable backrest
- Single-hook IV pole
- Removable side rails
- Anti-microbial mattress
- Set of 4 locking casters

## Specifications

- 500 lb. working capacity
- Sleeping Surface: 31" x 81"
- Height: 19" - 32" (adjustable)
- Stored: 18" x 31" x 3"
- Weight: 64 lbs.
- ADA compliant



info@VeriCorMed.com  
(866) 469-6019

For more information & pricing: [www.VeriCorMed.com](http://www.VeriCorMed.com)

**VERICOR** LLC  
Response-Friendly™