

Functional Needs Bed

Model: RB-FN400

The Functional Needs (FN) Bed is designed with patient surge capacity care needs in mind. Meeting FEMA guidelines and ADA compliant, this bed offers the patient a sense of safety and comfort. It has a wheelchair transferable height and a decontaminable mattress made to reuse over extended periods of time.



* IV pole and side rails options

Folded (Dual View)

Features

- Lightweight and portable
- 5 position head elevator
- 2 position foot elevator
- Anti-microbial mattress
- Spring and deck design
- Patient restraint belt
- Heavy duty locking legs

Specifications

- 450+ lbs. working capacity
- Sleeping Surface: 32" x 80"
- Height: 18"
- Folded: 34" x 42" x 8"
- Weight: 32 lbs.
- Meets FED STD 1633 for flammability



Options	Model #
w/ IV Pole	RB-FN400I
w/ IV Pole & Rails	RB-FN400ISR

[Go to Web Page >>](#)



info@VeriCorMed.com
(866) 469-6019

For more information & pricing: www.VeriCorMed.com

