42" Workstation w/ 3 Drawers, Table & Refrigerator

Model: MW-42R

VeriCor's 42" Workstation with Refrigerator, 3 Drawers & Table comes standard with ample storage, a built-in refrigerator and a detachable table. This single unit station becomes mobile in an instant. Pack in all your medical supplies and the workstation is ready for your next vaccination campaign.

Integrated Refrigerator







Locking Lid Stays



Swivel Wheels



Features

- Refrigerator without freezer
- Lockable refrigerator & case
- In-lid glove holders
- Easy lock wheels with a wide wheel base for easy and safe maneuverability
- Custom labeling included

Specifications

Dimensions (H x W x D)
48" x 43" x 25"

Drawers
3 (100 lbs. capacity each)

Fold-Out Table
35" x 43" (30" tall)

Weight
283 lbs.

Please provide color code when ordering:

COLOR (CODE)	Red (RD)	Yellow (YL)	Forest Green (FG)	Navy Blue (NB)	Black (BK)
Silver (SL)	Orange (OR)	Green (GR)	Violet Blue (BL)	Wine (WN)	Tan (TN)

Go to Web Page >>

info@VeriCorMed.com
(866) 469-6019

Made In the U.S.A.

SS-MW-42R-091026

For more information & pricing: **www.VeriCorMed.com**

VERICOR LLC
Response-Friendly®